




The Process of AFS

Stage 2. The radiation at 254nm is absorbed by any analyte mercury atoms and is emitted as fluorescence at 254nm also.




Element Specific
Selectivity




PSA

AFS - How do we capture the emitted photons? The importance of detector design – measurement at 90°

AFS



AAS



Emission v Absorption

Emission Technique
AFS/ICP

Absorption Technique
AAS/GF

Standby
(Background)



Analysis
(Signal)




Great sensitivity



PSA



Linearity



Emission techniques such as AFS

Several orders of magnitude linearity



PSA

Interferences

- Spectral (polyatomic - isobaric)
- Contamination and Carryover
(memory effects)

Spectral interferences

The Achilles heel of ICP-MS?
Consider Arsenic and Selenium

Table 1. Common ICP-MS interferences.

Polyatomic Species	Interfered Analyte
$^{40}\text{Ar}^{35}\text{Cl}$	^{75}As
ArAr	^{80}Se

Arsenic ^{75}As 100%
Selenium ^{74}Se 0.89%
 ^{76}Se 9.37%
 ^{77}Se 7.63%
 ^{78}Se 23.77%
 ^{80}Se 49.61%
 ^{82}Se 8.73%

Choose another isotope?

Use a high resolution MS?

Use He as a Plasma gas?

Use a collision cell?

There are no such
interferences for AFS



PSA

Contamination and Carryover (memory effects)

- Hg is notorious for causing carry over problems – Sticks to most surfaces – Care is required.
- A major problem is the presence of glassware within the instrumentation
 - Spray chamber, nebuliser, flow cells



PSA – AFS – Patented
Argon shield flow cell

